


MALDI Biotyper – Changing Microbiology




4. Freiburger MALDI-Meeting – 25 Juni 2019

Dr. Gerold Schwarz - Sales and Support - Bruker Microbiology & Diagnostics Bremen

Bruker Corporation and MALDI Biotyper



- 1960 Foundation of the **Bruker** company
- 1992 Introduction **first MALDI-TOF** mass spectrometer
- 2004 Launch of **MALDI Biotyper** as a research tool
- 2008 First MALDI Biotyper in **routine laboratories**
- 2009 Launch of the MALDI Biotyper **IVD system**
- 2013 **FDA Clearance** of the MALDI Biotyper CA system
- 2017/18 **AOAC**-approvals, **ISO 16140**-part 6 validation



1988 → MALDI technique by Karas M., Hillenkamp F. *Anal Chem.* 1988 Oct 15; 60 (20): 2299-2301



2002 → Nobel Prize in chemistry (for soft ionization techniques)

MALDI Biotyper

Microorganism Confirmation/Identification by MALDI-TOF MS: Workflow




Selection of colony of
unknown microorganism



Preparation onto
MALDI target plate

96 samples identified
in ~60 minutes



Identification by
pattern matching





Acquisition of MALDI-TOF
profile spectrum



MALDI Biotyper - Basics

Matrix **A**ssisted **L**aser **D**esorption / **I**onization
Time **O**f **F**light Mass Spectrometry



MALDI Biotyper

Sample preparation procedures



Direct transfer (for 90-95% of microorganisms)

- Direct transfer of biological material
- add MALDI Matrix
- Gram negative and most Gram positive microorganisms


Extended direct transfer (e.g. for yeast)

- Direct transfer of biological material
- Add drop of formic acid
- Dry at room temperature
- Add matrix

Extraction (for MSPs, better ID..)

- Enrich "targeted" molecules from isolated colonies

→ Special protocols for special organisms



Direct transfer using a tooth pick

**Robust pattern matching to
MBT Compass Library**

Updates on Libraries and Modules

MALDI Biotyper

Reference library



Bruker MBT Compass Library (version April 2019)

8468 MSPs

2950 different species

	MSP	Genus	Species
Gram -	3525	258	1266
Gram +	4070	210	1450
Yeast	806	45	210
Filamentous Fungi	67	27	43
Σ	8468	540	2969

Every Microorganism
could be of Relevance

Abiotrophia defective *Acetobacter* *aceti* *Acetobacter cerevisiae* *Acetobacter*
Achromobacter xylosoxidans *Acidaminococcus fermentans* *Acidaminococcus*
Acinetobacter bouvetii *Acinetobacter calcoaceticus* *Acinetobacter gementii* *Acinetobacter*
Acinetobacter schindleri *Acinetobacter sp* *Acinetobacter tandoii* *Acinetobacter*
Actinobacillus ureae *Actinobaculum* *massiliense* *Actinobaculum schaalii*
denticolens *Actinomyces europaeus* *Actinomyces funkei* *Actinomyces ge*
Actinomyces nasicola *Actinomyces neuii* *Actinomyces odontolyticus* *Actinomyces*
viscosus *Actinomyces weissii* *Adlercreutzia* *equolifaciens* *Advenella* *ince*
caviae *Aeromonas encheleia* *Aeromonas enteropelogenes* *Aeromonas eucr*
sobria *Aeromonas sp*[2] *Aeromonas veronii* *Afipia* *broomeae* *Afipia felis*
humatus *Agromyces italicus* *Agromyces lapidis* *Agromyces medolanus* *Ag*
contaminans *Alicyclobacillus cycloheptanicus* *Alicyclobacillus fastidiosus*
vulcanalis *Aliivibrio fischeri* *Alishewanella fetalis* *Alistipes finegoldii* *Alis*
coloradensis *Amycolatopsis fastidiosa* *Amycolatopsis japonica* *Amycolatop*
Anaerobiospirillum *succiniciproducens* *Anaerococcus* *hydrogenalis* *An*
migulanus *Aquicola* *tertiaricarbonis* *Arcanobacterium* *canis* *Arcanobac*
nitrofigilis *Arcobacter skirrowii* *Aromatoleum* *alkani* *Aromatoleum anae*
Aromatoleum toluvorans *Arsenicicoccus bolidensis* *Arsenicicoccus derm*
citreus *Arthrobacter creatinolyticus* *Arthrobacter crystallopoietes* *Arthrob*
Arthrobacter monumenti *Arthrobacter mysorens* *Arthrobacter nasiphocae* *A*
protophormiae *Arthrobacter psychrolactophilus* *Arthrobacter psychrophena*
tumbae *Arthrobacter uratoxydans* *Arthrobacter ureafaciens* *Arthrobacter*
Avibacterium gallinarum *Avibacterium volantium* *Azoarcus* *communis* *Azo*
arsenicus *Bacillus asahii* *Bacillus atrophaeus* *Bacillus azotoformans* *Bacillu*
Bacillus cohnii *Bacillus decolorationis* *Bacillus drementensis* *Bacillus endophy*
hemicellulosilyticus *Bacillus horikoshii* *Bacillus horneckiae* *Bacillus horti* *E*
Bacillus mannanyticus *Bacillus marisflavi* *Bacillus megaterium* *Bacillus m*
Bacillus pseudofirmus *Bacillus pseudomycoides* *Bacillus psychrosaccharo*
Bacillus subterraneus *Bacillus subtilis* *Bacillus thermoamylovorans* *Bacillu*
coagulans *Bacteroides coprocola* *Bacteroides coprophilus* *Bacteroides eg*
plebeius *Bacteroides pyogenes* *Bacteroides salyersiae* *Bacteroides sterco*
angulatum *Bifidobacterium animalis* *Bifidobacterium asteroides* *Bifidobacter*
longum *Bifidobacterium magnum* *Bifidobacterium merycicum* *Bifidobac*
Bifidobacterium thermophilum *Bifidobacterium* *Bifidobacterium* *Bifidobacterium*

MALDI Biotyper

Libraries new entries



ECCMID 2019:

- **Standard 9.0 8468 MSPs** (8326 IVD)
- **Mycobacteria 6.0** – (952 MSPs + IVD identical)
- **Filamentous Fungi 3.0** (new - 577 MSPs)
- SR library (IVD extension - BBFV)
- Third claim US-FDA (July 2017 - 424 bacteria)
- April 2018 US-FDA (candida auris added – 425 bacteria)
- Third parties (MALDI-Up, D-MASS, CDC-Microbenet, ..)

One workflow - Two calculations



- 2016

- 2017

Listeria monocytogenes species confirmation

MBT Subtyping Module

Listeria monocytogenes



- *Listeria monocytogenes* group members are closely related (*L. monocytogenes*, *L. ivanovii*, *L. innocua*, *L. welshimeri* and *L. seeligeri*)
- **Mainly *L. monocytogenes* is pathogenic for humans**
- *L. monocytogenes* is one of the most virulent foodborne pathogens, causing listeriosis
- The infection is most likely to sicken pregnant women and their newborns, adults aged 65 or older, and people with weakened immune systems

- **Rapid and accurate identification** of *Listeria* strains is essential for appropriate management and control of food safety



MBT Subtyping Module

Listeria monocytogenes

Result in MBT Compass: Report



- Listeria monocytogenes* typing to confirm species

→ **directly** from direct transfer (and extended direct transfer)

Sample Name	Sample ID	Organism (best match)	Score Value	Organism (second-best match)	Score Value
0_E1 (+++)(A)	15015734-1 (LM) TSA_ALOA 24 h Art EDT3 (standard)	<i>Listeria monocytogenes</i> typed as L. monocytogenes	2.30	Listeria monocytogenes	2.29
0_E2 (+++)(A)	15015734-1 (LM) TSA_ALOA 24 h Art EDT4 (standard)	<i>Listeria monocytogenes</i> typed as L. monocytogenes	2.31	Listeria monocytogenes	2.30

ID performed on DT => ID > 2.0 => Subtyping performed

NO additional steps required

MBT Subtyping Module

*bla*_{KPC} plasmid detection in *K. pneumoniae*

Result in MBT Compass: Report



Sample Name	Sample ID	Organism (best match)	Score Value	Organism (second-best match)	Score Value
<u>0_A7</u> (+++)(A)	KP-M-93 chrom (standard)	Klebsiella pneumoniae presumptive <u>KPC</u>	<u>2.47</u>	Klebsiella pneumoniae	<u>2.38</u>
<u>0_A8</u> (+++)(A)	KP-M-94 chrom (standard)	Klebsiella pneumoniae presumptive <u>KPC</u>	<u>2.49</u>	Klebsiella pneumoniae	<u>2.40</u>
<u>0_A9</u> (+++)(A)	KP-M-94 chrom (standard)	Klebsiella pneumoniae presumptive <u>KPC</u>	<u>2.48</u>	Klebsiella pneumoniae	<u>2.37</u>

Subtyping Result	Comment
KPC	classification as KPC on the basis of specific marker peaks (11109)

ID performed on DT => ID > 2.0 => Subtyping performed

NO additional steps required


MBT Subtyping Module

*bla*_{KPC} plasmid detection in *K. pneumoniae*



*bla*_{KPC} plasmid peak detection at m/z 11,109

- Good detection rate
- Very specific peak



Klebsiella pneumoniae:


- BLUE: KPC with *bla*_{KPC} related peak
- RED: ESBL

Automation

MBT Pilot and Galaxy

MALDI Biotyper

Paperless and traceable workflow



MALDI Biotyper

Paperless and traceable workflow

MBT Satellite software



For convenient sample prep on the work benches:

- Software module for project creation and set-up of corresponding working lists



- Also optionally available pre-installed on a tablet-PC (Display # 9.7")

MBT Pilot


Traceability, safety and confidence



- Complements the MBT Satellite software to a **barcoded and paperless workflow**
- **Micro-projection technology**
 - For optically guided preparation
 - No laser reflection
 - Patent pending
- **Cross hairs** to indicate next position
- Preparation spot **not covered**: no glare
- Process is fully **traceable**
- **Ergonomic**
- Tool for improved workflow and quality control

MBT Pilot


Guided Target Preparation



1. Scan barcode of MALDI target




2. Insert target into MBT Pilot target holder, the next free position is indicated on the target



3. Apply sample material at indicated position




4. Scan sample barcode



5. Prepare additional samples by repeating steps 3 and 4

MALDI Biotyper

Subsequent sample preparation steps



MBT Galaxy

For automated application of matrix and formic acid





- Frees time for laboratory personnel
- Ensures highest preparation quality
- Improves traceability
- Less waste

- Frees laboratory personnel from pipetting
- Onboard barcode reading automatically loads the work list from MBT server
- Automated formic acid / matrix preparation
- Contact free droplet application
 - ✓ eliminates carry-over
 - ✓ no pipette tips → less waste
- Ensures highest preparation quality by control of:
 - Previous matrix preparation on the spot
 - Droplet volume
 - Drying atmosphere
- Improves traceability in a paperless workflow

MBT Pilot and MBT Galaxy

Combined workflow




1. Press Start

→ Insert target



2. Select preparation

→ Target identification
→ Spotlist incl. BTS
submitted by MBT Server



3. Target preparation

MALDI Biotyper

Combined workflows - other modules



- Fungi Library 3.0 and MBT filamentous Fungi Module
- Mycobacteria 6.0 and MBT Mycobacteria Module
- Subtyping Module
- Food Report Option (AOAC, MicroVal, Confirmation)
- Sepsityper Module (direct blood culture analysis)
- STAR BL Module (Selective Testing of Antibiotic Resistance)
- IVD- mode, RUO- mode, GP (industrial microbiology),

Automation

Colony picker, robotic..


MALDI Biotyper

Micronaut ASTroid



MICRONAUT ASTroid

Workflow MICRONAUT /
MALDI-TOF fully
automated



Inoculation of MICRONAUT MIC plates and MALDI-TOF targets from a uniform bacteria suspension

MALDI Biotyper

BD Kiestra Lab Automation



IDENTIFICATION | MALDITOF



Automatic Colony Picking

Transferring colonies to a MaldiTOF target plate is both labour intensive and error sensitive. It is now possible to automate this vital process with the MaldiTOF, which helps you to improve quality and efficiency of the identification step. The automatic preparation of the target plate will ensure that the right amount of sample material and the right amount of matrix and acid fluids will be added for direct and extended-direct transfer.



Automatic Suspension preparation

The system will automatically prepare a suspension fluid for AST with a desirable McFarland concentration. An integrated density meter will ensure the right concentration is prepared. The prepared AST suspension tubes can be placed in the InoQuA for automatic inoculation and spreading. This will result in a perfect and fast prepared confluent spread.

MALDI Biotyper MALDI Colonyst..




MALDI Colonyst[®]




MALDI Biotyper


MALDI Colonyst..




MALDI target and Petri dish inserted for transfer of colonies




Pre-deposition of formic acid using liquid deposition needle




Formic acid deposited only on preselected positions



Each selected colony precisely picked...



...and smoothly deposited into formic acid droplet



After drying of all colonies, MALDI matrix is deposited

with feeder



MALDI Biotyper

Latest News



Bruker Launches MALDI Biotyper Sirius at ASM Microbe Conference

... new MALDI Biotyper *sirius*[™] system with additional negative-ion mode to support MALDI research and RUO validation studies, e.g. for colistin-resistance testing

US introduction of Micronaut[™] products for true minimum inhibitory concentration (MIC) antibiotic susceptibility testing (AST) in veterinary medicine

US introduction of RUO versions of MALDI Biotyper-based rapid Selective-Testing of Antibiotic Resistance (MBT-STAR[™]) assays for validation studies on carbapenem resistance and cephalosporine resistance






The new RUO colistin-resistance assay uses lipid analysis in negative-ion mode, and has been developed at Imperial College London, UK

Larrouy-Maumus et al., presented at ECCMID 2019, (submitted for publication).

Article | [OPEN](#) | Published: 15 November 2018

Rapid detection of colistin resistance in *Acinetobacter baumannii* using MALDI-TOF-based lipidomics on intact bacteria

Laurent Dortet , Anais Potron, Rémy A. Bonnin, Patrick Plesiat, Thierry Naas, Alain Filloux & Gerald Larrouy-Maumus 

Scientific Reports **8**, Article number: 16910 (2018) | [Download Citation](#) 

MALDI Biotyper

Official Method of Analysis by AOAC International



- The MBT has been certified according to the

Official Method of Analysis program
(OMA) of the AOAC International

for the **confirmation and identification**
of:

- ✓ *Salmonella* spp.
Cronobacter spp.
Campylobacter spp.
and other gram-negative bacteria
- ✓ *Listeria* spp.
Listeria monocytogenes
and other gram-positive bacteria



First Action AOAC Official MethodSM
2017.09

First Action AOAC Official MethodSM
2017.10

MALDI Biotyper

ISO 16140-Part 6 Certification by MicroVal



- The MBT is the very **FIRST and ONLY** confirmation method certified by MicroVal according to the NEW **ISO 16140-part 6 standard**



for the **confirmation of:**

- ✓ ***Cronobacter spp.***
- ✓ ***Salmonella spp.***
- ✓ ***Campylobacter spp.***
- ✓ ***Listeria spp. and Listeria monocytogenes***

Certificate N° 2017LR72

Certificate N° 2017LR73

Certificate N° 2017LR74

Certificate N° 2017LR75

Thank you for your attention!

Dr. Gerold Schwarz
Applications Support Team Bremen

gerold.schwarz@bruker.com
biotyper.appl.support.EMEA@bruker.com



Bild: Bruker Daltonik GmbH Bremen